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Cat: G1281

Size: 4×50mL/4×100mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent	4×50mL	4×100mL	Storage
Reagent (A): PAS Oxidant	50mL	100mL	2-8°C, avoid light
Reagent (B): Schiff Reagent	50mL	100mL	2-8°C, avoid light
Reagent (C): Mayer Hematoxylin Solution	50mL	100mL	RT, avoid light
Reagent (D): Acidic Differentiation Solution	50mL	100mL	RT

Introduction

Glycogen staining is one of the conventional staining methods in pathology. McManus first used PAS technology to display mucin in 1946. This method is often used to display glycogen and other polysaccharides. The staining solution can not only display glycogen, but also show neutral mucilaginous substances and some acidic substances, as well as cartilage, pituitary, mould, fungus, pigment, amyloid substance, basement membrane, etc.

The oxidant can oxidize 1,2-glycol group in sugars and related substances to make it into dialdehyde. Aldehyde and Schiff Reagent can combine to form a purplish-red fuchsin complex. As oxidants can also oxidize other substances in cells, we should pay attention to the concentration and oxidation time of oxidants when using, so that the oxidation control can oxidize glycol group to aldehyde group without peroxide, which is a key step.

This product has some advantages: the unique formula technology of Solarbio greatly enhances the dyeing effect; stable performance and strong specificity; simple operation and short time.

Self Provided Materials

10% Neutral Buffered Formalin Fixative, Alcohol Acetate Formalin(AAF) Mixed Fixative(70%) or Carnoy Fixative, distilled water, ethanol.

Protocol(for reference only)

1. Carnoy or AAF fixative is recommended for glycogen staining, and 10% formalin fixative is recommended for glycoprotein and glycolipid, and routine dehydration and embedding is recommended. (See Note 2)
2. Conventionally dehydration and embedding. For paraffin section, dewax to distilled water; For frozen section, directly into distilled water to restore room temperature.
3. Place in Oxidant for 5-8min at room temperature, generally not more than 10min.
4. Wash with distilled water twice, each time for 1min.
5. Soak the section in Schiff Reagent in a dark place at room temperature, and stain for 10-20min.
6. Rinse with distilled water for 1-2min.
7. Stain with Mayer Hematoxylin Solution for 1-2min.
8. (Optional) Differentiate by Acidic Differentiation Solution for 2-5s.
9. Wash with tap water for 10-15 min or rinse in 1×PBS for 2-3min to blue.
10. Starting from 75% ethanol, dehydrate with 85%, 95%, and 100% ethanol step by step, with each step lasting 1-2 minutes. Transparent by xylene twice, each time for 1 minute, and then seal with resinene.
11. View under the microscope and photograph.

Result

PAS Reaction Positive Substance	Red or Purplish Red
Nucleus	Blue
Cytoplasm	Red in different degrees

Note: the color depth depends largely on the length of time that the sample has been in the Oxidant and Schiff Reagent.

Negative Control (optional)

1. Take 1gm of amylase and dissolve it in 100mL of PBS (pH5.3), treat it for 30-60 min, and add into Oxidant together with other sections. The result should be negative.

2. (alternative) take the saliva(after filtration) and treat it for 30-60 min, then add into Oxidant together with other slices. The result should be negative.
3. (alternative) if the control section use its own sample, directly place the section into Schiff Reagent without Oxidant. The result should be negative.

Note

1. Section dewaxing should be as clean as possible, otherwise it will affect the dyeing effect.
2. Glycogen is not tightly bound to tissues, and usually does not cross-link with aldehydes. It is recommended to use alcohol-based fixative with no or less water for fixation. Glycoproteins or glycolipids can be firmly bound to tissues or cross-linked and fixed with aldehydes, so conventional aldehyde fixative can be used.
3. The oxidation time of PAS Oxidant should not be too long, and the best temperature is 18-22 °C.
4. The Oxidant and Schiff Reagent should be kept in an airtight place at 4 °C to avoid too much sunlight and air during use. Before use, it is better to take it out 30 min in advance to restore room temperature, and use it in dark.
5. Acid Differentiation Solution should be replaced frequently, and its differentiation time should be determined according to the thickness of section, the type of tissue and the old and new of acid differentiation solution. In addition, the washing time of tap water after differentiation should be enough.
6. The time of action in Oxidant and Schiff Reagent is very important, which depends on the thickness of section and the type of tissue.
7. This product is often used for conventional tissue section staining. For fungi, cells and extremely thin sections, it is recommended to purchase Glycogen Periodic Acid Schiff (PAS) Stain Kit (For Fungus), because the concentration of Oxidant and Hematoxylin Solution is lower, so it is not suitable for over staining.
8. The staining time of frozen section should be as short as possible.
9. For your safety and health, please wear experimental clothes and disposable gloves.