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Alkaline Phosphatase Stain Kit (Modified Gomori Ca-CoS Method)

Cat: G1481

Size: 3×50mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent	3×50mL	Storage
Reagent(A): ALP Incubation Solution	50mL	2-8°C, avoid light
Reagent(B): Co Solution	50mL	RT, avoid light
Reagent(C): Vulcanizing Solution	2×1mL	RT, avoid light
Dilute Reagent(C) 50 times with distilled water or deionized water to obtain the Vulcanization Working Solution, which is prepared and used immediately.		
Reagent(D): ALP Contrast Solution	10mL	2-8°C, avoid light

Introduction

Alkaline phosphatase (ALP or AKP) is a kind of phosphatase widely distributed in mammalian tissues. The optimum pH for its activity is 9.2-9.8. The phosphatase mainly exists in the active sites of substance exchange (cell membrane), such as brush-like margin of intestinal epithelium and proximal convoluted tubule of kidney, stationary cilia of epididymis epithelium, capillary bile duct membrane of liver, and endothelium of arterioles and capillary arteries. It also exists in endoplasmic reticulum, Golgi complex, ingestive vesicles, lysosomes of intestinal epithelium, neutrophils, cell membranes and ingestive vesicles of smooth muscle.

This kit uses metal precipitation method to display the alkaline phosphatase activity. Alkaline phosphatase hydrolyzes to release phosphoric acid using naturally occurring sodium beta-glycerophosphate as the substrate. Then the phosphoric acid is precipitated into calcium phosphate by calcium ion immediately, replaced with cobalt phosphate by Co Solution again, and finally replaced with black precipitation by Vulcanizing Solution.

Self Provided Materials

Distilled water, Incubator, Water bath, 1% calcium chloride aqueous solution

Protocol (for reference only)

For paraffin section(See Note 3)

1. Dewax paraffin sections and rehydrate in graded alcohol to distilled water.
2. Cover the section with Reagent(A): ALP Incubation Solution and incubate at 37°C for 2-12 h.
3. Treat with 1% calcium chloride three times, each time for 1 minute. After each treatment, pour out excess liquid and do not rinse vigorously. (See Note 4)
4. Cover the slices with Reagent (B): Co solution and incubate at 37 °C for 5 min.
5. Rinse with running water for 5min and place it in distilled water.
6. Dilute reagent (C):Vulcanizing Solution 50 times with distilled water to prepare the Vulcanization Working Solution, which is ready to use. (See Note 2)
7. Cover the section with Vulcanizing Working Solution and incubate for 1-2min till it turn black.
8. Rinse with distilled water for 3 times, each time for 3 min.
9. (Optional) Re-dyeing with Nuclear Fast Red for 3min and wash with distilled water for 1min.
10. Absolute ethanol dehydration, transparent by xylene, seal with resinene.

For frozen section(See Note 3)

1. Fix the frozen section by acetone for 2-5min .
2. Cover the section with Reagent(A): ALP Incubation Solution and incubate at 37°C for 45-75 h.
3. Treat with 1% calcium chloride three times, each time for 1 minute. After each treatment, pour out excess liquid and do not rinse vigorously. (See Note 4)
4. Cover the slices with Reagent (B): Co Solution and incubate at 37 °C for 5 min.
5. Rinse with running water for 5min and place it in distilled water.
6. Dilute reagent (C):Vulcanizing Solution 50 times with distilled water to prepare the Vulcanization Working Solution, which is ready to use. (See Note 2)
7. Cover the section with Vulcanizing Working Solution and incubate for 1-2min till it turn black.
8. Rinse with distilled water for 3 times, each time for 3 min.

9. (Optional) Re-dying with Nuclear Fast Red for 3min and wash with distilled water for 1min.
10. Seal with Glycerol gelatin or usually dehydrate, transparent by xylene, seal with resinene.

Result

Positive Site	Black Precipitation
Nucleus	Red(By Nuclear Fast Red)

Negative Control(optional)

1. Reagent negative: Reagent (D) is a substrate free incubation solution. Take adjacent slices and add Reagent (D): ALP Contrast Solution substitute Reagent(A):ALP Incubation Solution for incubation. The staining result is negative for the reagent.
2. Negative enzyme activity: Before staining the sections, alkaline phosphatase was inactivated in iodine solution and 5% sodium thiosulfate solution for 3 min each. After thorough washing with water, incubation and other steps were carried out. The staining result was negative for enzyme activity.

Note

1. Reagent (A): ALP Incubation Solution is a sterile solution. It is recommended to receive the reagent kit and store it in small portions on the ultra clean table. Once opened, it should be used immediately.
2. Reagent(C): Vulcanizing Solution has a special odor and is corrosive. Use and dilution operations should be carried out carefully in a fume hood.
3. Alkaline phosphatase is relatively fragile, and samples should be processed immediately after sampling. It is recommended to use a refrigerator at 2-8 °C for a fixed period of time, which should not exceed 24 hours. The recommended embedding temperature should not exceed 60 °C and the time should not exceed 3 hours to avoid inactivation leading to no positive results.
4. Due to the formation of calcium phosphate that is slightly soluble in water, washing with distilled water may wash away the calcium phosphate, resulting in subsequent staining without positive or weak positive results. Gentle cleaning with 1% calcium chloride will replenish calcium ions to ensure the stability of calcium phosphate.
5. For your safety and health, please wear laboratory clothes and disposable gloves when operating.