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- [1] Dan Xiao,Yue Zhang,Rui Wang,et al. Emodin alleviates cardiac fibrosis by suppressing activation of cardiac fibroblasts via upregulating metastasis associated protein 3. Acta Pharmaceutica Sinica B. April 2019. (IF 5.808)
- [2] Na Li,Lin Zhou,Weilong Xie,et al. Alkaline phosphatase enzyme-induced biomineralization of chitosan scaffolds with enhanced osteogenesis for bone tissue engineering. Chemical Engineering Journal. September 2019;371:618-630. (IF 8.355)
- [3] Ke Xue,Jun Zhang,Cong Li,et al. The role and mechanism of transforming growth factor beta 3 in human myocardial infarction-induced myocardial fibrosis. Journal of Cellular and Molecular Medicine. April 2019. (IF 4.658)

Modified Masson's Trichrome Stain Kit

Cat: G1346

Size: 8×50mL/8×100mL

Storage: 2-8°C, avoid light, valid for 1 year.

Kit Components

Reagent	8×50mL	8×100mL	Storage
Reagent(A): Mordant Solution(3×)	50mL	100mL	RT, avoid light
Reagent(B): Celestite Blue Solution	50mL	100mL	2-8°C, avoid light
Reagent(C): Mayer Hematoxylin Solution	50mL	100mL	2-8°C, avoid light
Reagent(D): Acid Differentiation	50mL	100mL	RT
Reagent(E): Ponceau-Acid Fuchsin Solution	50mL	100mL	RT, avoid light
Reagent(F): Phosphomolybic Acid Solution	50mL	100mL	RT, avoid light
Reagent(G): Aniline Blue Solution	50mL	100mL	RT, avoid light
Reagent(H): Weak Acid Solution(3×)	50mL	100mL	RT

Introduction

Narrowly speaking, connective tissue refers to the three types of fibers it contains: collagen fibers, reticular fibers, and elastic fibers. Collagen fibers are the most widely distributed and abundant type of fiber. Masson trichrome staining, also known as Masson staining, is the most classic method in connective tissue staining and an authoritative and classic technical method for collagen fiber staining. The so-called trichrome staining usually refers to staining the cell nucleus and selectively displaying collagen and muscle fibers. The dyeing principle of this method is related to the size of anionic dye molecules and the permeability of tissues: the size of molecules is reflected by their molecular weight. Small molecular weights are easy to penetrate tissues with dense structures and low permeability, while large molecular weights can only enter tissues with loose structures and high permeability. However, the molecular weight of light green or aniline blue is large, so after Masson staining, muscle fibers appear red, while collagen fibers appear green or blue, mainly used to distinguish between collagen fibers and muscle fibers.

The difference between the Modified Masson's Trichrome Stain Kit and the conventional Masson Trichrome Stain Kit is that the nucleus of the cell is light stained with Celestite Blue Hematoxylin Solution. It has the following characteristics: short differentiation time, clear and bright colors; wide application range, suitable for paraffin section, frozen section and other staining of tissues; long preservation time of stained sections and not easy to fade. Modified Masson's Trichrome Stain Kit shows that collagen fiber is blue, muscle fiber, cytoplasm, cellulose, keratin and erythrocytes are red, and nucleus is blue, which was mainly used to distinguish collagen fiber and muscle fiber.

Self Provided Materials

G2161-Neutral Buffered Formalin Fixative, 10% or P1110-Paraformaldehyde, 4%, Bouin fixative, distilled water, series ethanol, xylene, or G2150 environmental protection tissue transparent dewaxing solution, dye tank

Protocol(for reference only)

Before staining, prepare the working solutions according to the required amount as follows: Mix Reagent (A): Mordant Solution (3×) with distilled water at a ratio of 1:2 to obtain Mordant Solution (1×). Mix Reagent (H): Weak Acid Solution (3×) with distilled water at a ratio of 1:2 to obtain Weak Acid Working Solution.

For paraffin sections

1. Fix tissues in 10% neutral formalin or Paraformaldehyde, 4%, then dehydrate and embed.
2. Cut the section in 3-6μm thick, dewax to distilled water.
3. Incubate the section in Mordant Solution(1×) for 8-12h at RT or incubate the section in 57-60°C incubator for 1h, then wash with distilled water for three times, each for 3min. (See Note 2)
4. Drop Celestite Blue Solution onto the section and stain for 2-3min. Slightly wash with distilled water twice, each for 10-15s.
5. Drop Mayer Hematoxylin Solution onto the section and stain for 2-3min. Slightly wash with distilled water twice, each for 10-15s.
6. Differentiate by Acid Differentiation for several seconds and wash with distilled water to stop differentiation.

Rinse in running water for 10min(Return to Blue). (See Note 2)

7. Drop Ponceau-Acid Fuchsin Solution onto the section and stain for 5-10min. Slightly wash with distilled water twice,each for 10-15s .
8. Treat with Phosphmolybic Acid Solution for 5-10min.
9. Discard the remaining dye solution without washing and stain with Aniline Blue Solution for 3-5min.
10. Wash the section with Weak Acid Working Solution, and rinse in Weak Acid Working Solution for 2mins.
11. Dehydrate quickly in 95% ethanol. Dehydrate in absolute ethanol twice for each time 5-10s.(See Note 6)
12. Transparent by xylene twice for each time 1-2mins. Seal with resinene.

For frozen sections

1. Fix tissues in 10% neutral formalin or Paraformaldehyde,4%, dehydrate by gradient sucrose and embed in OCT to prepare frozen sections.
2. Cut the section in 8-15 μ m thick. Before staining, place the slices in distilled water and restore to room temperature for 3min.
3. The staining steps are the same as steps 3-6 for paraffin sections.
4. Drop Ponceau-Acid Fuchsin Solution onto the section and stain for 5-7min. Slightly wash with distilled water twice,each for 10-15s .
5. Treat with Phosphmolybic Acid Solution for 2-3min.
6. Discard the remaining dye solution without washing and stain with Aniline Blue Solution(diluted 3-4 times by Weak Acid Working Solution)for 1-3min. (See Note 5)
7. Wash the section with Weak Acid Working Solution, and rinse in Weak Acid Working Solution for 2mins.
8. Dehydrate quickly in 95% ethanol. Dehydrate in absolute ethanol twice for each time 5-10s.(See Note 6)
9. Transparent by xylene twice for each time 1-2mins. Seal with resinene.

Result

Collagen Fiber	Blue
Muscle Fiber, Cytoplasm, Cellulose, Keratin and Erythrocytes	Red
Nucleus	Bluish Brown

Note

1. Section dewaxing should be as clean as possible.
2. The Mordant Solution and Acid Differentiation of this kit are optimized formulas, which are colorless solvents without dangerous chemicals. They are quite different from the properties of classical picric acid mordant and acid ethanol differentiation solution, witch belongs to normal phenomena. The above two reagents are optimized substitutes for the original reagent. If you prefer the classic formula, you can replace the reagents used in the corresponding steps by yourself.
3. The differentiation time of Acid Differentiation should be determined by the thickness of section, the type of tissue and the old and new.
4. Phosphmolybic Acid Solution has two effects. One is that Phosphmolybic Acid Solution can make the collagen fiber dyed red differentiate into colorless or light red, while the muscle fiber cellulose is still bright red; the other one is that it can also play a role of mordant dyeing on the collagen fiber, so that the collagen fiber is easy to combine with Aniline Blue Solution of macromolecular dye.
5. Compared to paraffin sections, frozen sections are prone to staining collagen fibers with aniline blue. Therefore, it is recommended to dilute the Aniline Blue Solution with a Weak Acid solution 3-4 times before staining. For skin and other easily stained tissues, the staining time can be shortened again.
6. If the tissue shows red or blue decolorization after treatment with 95% ethanol and absolute ethanol, after cleaning the tissue with the Weak Acid Working Solution, you can spin off the excess liquid, let it air-dry directly, and then proceed with the subsequent xylene clearing and mounting steps.
7. Weak Acid Solution can make the color more clear and bright. If the use amount is large, you can buy G2940-Weak Acid Solution instead.
8. For your safety and health, please wear experimental clothes and disposable gloves.