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**Cat:** G3061

**Size:** 20mL/100mL

**Storage:** -20°C, avoid light, valid for 6 months.

### Kit Components

Reagent	20mL	100mL	Storage
Reagent(A1): X-Gluc Powder	14mg	2×35mg	-20°C, avoid light
Reagent(A2): X-Gluc Solvent	0.4mL	2×1mL	RT, avoid light
Reagent(B): GUS Buffer	20mL	2×50mL	2-8°C, avoid light

### Introduction

X-Gluc(X-GlcA) has a molecular weight of 521.8 and CAS number of 18656-96-7. It is a substrate for the detection of GUS gene in *Bacillus coli* and can rapidly detect GUS gene fusion markers in plants. Under suitable reaction conditions,  $\beta$ -glucosidase (GUS) can hydrolyze X-Gluc into blue substance, which is not dissolved in Indigo substance in transgenic nuclear tissue. The sites or sites with GUS activity show blue or blue spots, which can be observed by naked eye or microscope. The GUS Staining Kit is mainly used for GUS gene expression analysis of transgenic plants.

### Protocol(for reference only)

1. Prepare X-Gluc Solution(50×): For the 20mL specification, take 0.4mL of X-Gluc Solvent and add it into a 14mg X-Gluc Powder bottle to dissolve completely, which is the X-Gluc Solution(50×). For the 100mL specification, take 1 tube of 1mL X-Gluc Solvent and add it into a 35mg X-Gluc Powder bottle to dissolve completely, which is the X-Gluc Solution(50×). Then divide it into small size packages and store them at -20 °C or -80 °C. (See Note 2)
2. Prepare GUS Solution: take appropriate X-Gluc Solution(50×) and GUS Buffer, mix fully as the ratio of 1:50 to form GUS Solution. Such as mixing 200ul X-Gluc Solution(50×) and 10mL GUS Buffer to form 10mL GUS Solution. The solution is best to prepare before use and can store for 3 days at 4°C.
3. Take appropriate leaves and other tissues to be dyed and add appropriate GUS Solution to completely submerge the tissue.
4. Incubate at 37 °C for 1-24h. With the increase of incubation time, blue gradually appears. When the expression level is high, the sites or sites of GUS activity show blue or blue spots.
5. Remove the chlorophyll of the sample with 75% ethanol, and generally soak the sample in ethanol for 1-3h until the negative control shows white color. If necessary, repeat the decolorization step to completely remove chlorophyll. The samples are stored in ethanol and can be observed by naked eye or ordinary optical microscope. Blue on white background is GUS expression site.

### Note

1. The prepared GUS Solution can be stored for half a month at 4 °C in dark.
2. X-Gluc Solution(50×) should avoid repeated freezing and thawing, otherwise the dyeing efficiency will decrease.
3. Due to the tissue specificity and other reasons, the blue color reaction may not be completely consistent, so pay attention to explore the specific experimental conditions. The roots, flowers and leaves of *Arabidopsis* and the roots of tobacco seedlings can be dyed directly without any pretreatment. But the stems and leaves of plants like tobacco and potato must be cut into slices (1-3mm) before dyeing. When handling with large tissues and samples, vacuum infiltration method can be used to help substrate penetrate into cells. GUS Stain Kit (G3060) is recommended.
4. For your safety and health, please wear lab coats and disposable gloves when operating.