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规格：100mL/500mL

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- [1] Ke Xue,Jun Zhang,Cong Li,et al. The role and mechanism of transforming growth factor beta 3 in humanmyocardial infarction-induced myocardial fibrosis. Journal of Cellular and Molecular Medicine. April 2019.(IF 4.658)
- [2] Ying Shen,Aiping Qin. Regulation of Embryonic Signal on Talin1 in Mouse Endometrium. ReproductiveSciences. December 2018.(IF 2.548)
- [3] Danping Hong,Jiongyan Ding,Ouyang Li,et al. Human-induced pluripotent stem cell-derived macrophagesand their immunological function in response to tuberculosis infection. Stem Cell Research & Therapy.February 2018.(IF 4.963)

Giemsa Stain Solution (Working Suit)

Cat: G1010

Size: 100mL/500mL

Storage: RT, valid for at least 3 years.

Components

Reagent	100mL	500mL	Storage
Giemsa Stock Solution	10mL	50mL	RT
Giemsa Buffer	100mL	500mL	RT
Before use, mix Giemsa Stock Solution with Giemsa Buffer as the ratio of 1:9 to prepare Giemsa Working Solution, it is not advised to prepare in advance.			

Introduction

Giemsa Stain Solution is composed of azure and eosin. The chemical properties of various cell components are different, and the affinity for various dyes is also different. Eosinophilic particles are basic proteins, which are combined with eosin and dyed pink; eosinophilic particles, such as nucleoprotein or lymphocyte cytoplasm, are acid substances, which are combined with basic dye methylene blue or azure and dyed purple blue; neutral particles are in isoelectric state and can be combined with eosin and methylene blue, which are dyed light purple. All kinds of components are distinguished because of their own characteristics and the combination of different substances in Giemsa Stain Solution to present different colors.

The Giemsa Stain Solution prepared by our company is made of imported Giemsa dye. It can dye the nucleus to purplish red or blue purple, and the cytoplasm to pink, presenting a clear cell dyeing image under the light microscope. It is mainly used to show the difference of the shape and size of various blood cells in blood smear, with good staining effect, strong staining power and clear staining.

Protocol(for reference only)

1. Mix 1 part of Giemsa Stock Solution and 9 part of Giemsa Buffer to form Giemsa Working Solution. It is ready to use.
2. Prepared blood smear which focused on uniform cell distribution by routine method. Thoroughly dry blood or bone marrow smears. Fix smears in absolute methanol for 2-3 min.
3. Stain smears in Giemsa Working Solution for 15 min at room temperature.
4. Rinse slowly from one side of the section with distilled water or PB buffer with pH6.7-6.8 (be careful not to pour out the dye or rinse the membrane directly), dry and view under the microscope.

Note

1. For your safety and health, please wear experimental clothes and disposable gloves.
2. If the dyeing is too deep or too light, please adjust the dyeing time or the concentration of Giemsa Working Solution by yourself.
3. A small amount of precipitate may appear after diluting the Gimesa Stain, which will not affect the use; wash with water after dyeing.
4. Gimesa Stain is prepared by conventional method and filtered by filter paper. Please do not touch water when using. Otherwise, it will fail if the time is long.
5. PH has influence on cell staining. Slides for dyeing must be clean and free of acid and alkali pollution so as not to affect the dyeing results.

Related Products

G1011 Giemsa Stain Dilution Buffer
G1015 Giemsa Stain Solution(10*Stock Solution)
G1020 Wright-Giemsa Stain Solution
G1021 Wright-Giemsa Stain Kit
G1040 Wright Stain Solution
G4640 Fast Giemsa Stain Kit