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Cat : G1142

Size: 2×50mL/2×100mL/2×500mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	2×50mL	2×100mL	2×500mL	Storage
Reagent A: Alcoholic Hematoxylin Solution	50mL	100mL	500mL	RT, avoid light
Reagent B: Acid Ferric Chloride Solution	50mL	100mL	500mL	RT, avoid light
Mix equal parts of Reagent A and B to form Weigert's Hematoxylin Working Solution before use. This working solution is stable for about 24 hours.				

Introduction

Hematoxylin is one of the most commonly used dyes in histochemistry and immunohistochemistry. It can be widely used for nuclear staining of tissue sections or cultured cells. VG(Van Gieson)Staining is the first traditional fine staining method used to distinguish collagen fiber and muscle fiber, which is double stained. After staining, collagen fiber is bright red, muscle fiber, glial cell cytoplasm and red blood cell are yellow. Weigert's Hematoxylin Stain Kit is mainly used in the technology of acid treatment after nuclear staining, which can be used in combination with VG staining. After staining, the cell nucleus is blue or brownish blue.

Protocol(for reference only)

1. Dewax to distilled water.
2. Stain with Weigert's Hematoxylin Working Solution for 10-20 min. (Can be optimized and adjusted appropriately based on the pre experimental results). Rinse with distilled water to remove excess dye solution.
3. Treat with 0.5%-1% hydrochloric acid ethanol differentiation solution for 5-10 seconds, the tissue can be seen to turn slightly red with the naked eye. Rinse with distilled water to remove the differentiation solution.
4. Blue in tap water for 5-10 min.
5. Rinse with distilled water for 1 min.
6. (Optional) Stain with counterstain solution (modified VG staining solution) for 0.5-2 min.
7. (Optional) Quickly wash with distilled water for 5-10 seconds to remove excess dye solution.
8. Starting from 75% ethanol, rapidly dehydrate step by step, with 3-5 seconds for each. Transparent by xylene twice, each time for 1 min. Seal with neutral gum, and view and photograph under the microscope.

Result

Nucleus	Blue to Black
Collagen	Red
Muscle Fibers, Cytoplasm	Yellow

Note

1. When using this reagent for the first time, it is recommended to take 1-2 samples for pre-test. When there are a large number of samples, the dyeing rack and cylinder can be used for easy operation.
2. For your safety and health, please wear laboratory clothes and disposable gloves.