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Cat: G1285

Size: 6×50mL/6×100mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent	6×50mL	6×100mL	Storage
Reagent (A):Alcian Blue Staining Solution	50mL	100mL	2-8°C, avoid light
Reagent (B):PAS Oxidant	50mL	100mL	2-8°C, avoid light
Reagent (C):Schiff Reagent	50mL	100mL	2-8°C, avoid light
Reagent (D):Hematoxylin Solution	50mL	100mL	RT, avoid light
Reagent (E):Acidic Differentiation Solution	50mL	100mL	RT
Reagent (F):Scott Bluing Solution	50mL	100mL	RT

Introduction

Some cells, such as the Goblet cell of the gastrointestinal tract, can secrete thick secretions, which are mainly composed of Glycosaminoglycan and glycoprotein, and are also collectively referred to as mucus. Mucus can be divided into neutral and acidic mucus. In normal gastrointestinal mucus, the surface epithelium of the gastric mucosa, pyloric glands, and duodenal glands are mainly composed of neutral mucus substances; The Goblet cell and intestinal glands of the small intestine and large intestine mainly secrete acidic mucus. The combination of Alcian Blue (AB) and colorless fuchsin (PAS) staining can be used to distinguish neutral and acidic mucus, as well as their ratio. Some cancers can also be distinguished by AB-PAS staining, such as intestinal gastric cancer secreting acidic mucus in blue, and gastric cancer secreting neutral mucus in red.

This staining method is often used for the analysis of metaplasia degree of gastrointestinal lesions, and also for the identification of glandular secretions and Goblet cell. Our company integrates imported raw materials and optimized formulas into a ready-to-use dyeing kit, which can meet various dyeing needs for different purposes.

Self Provided Materials

Distilled water, Series of ethanol

Protocol (for reference only)

1. Dewax to distilled water, then rinse in distilled water for 2min.
2. Stain with Alcian Blue Staining Solution for 10-20min. Rinse in distilled water three times for each time 1-2min.
3. Treat by Oxidant for 5-8min. Rinse in distilled water twice, each time for 1-2min.(see note 3)
4. Soak in Schiff Reagent and stain for 10-20min. Then wash with distilled water twice.(see note 4)
5. (optional) Stain with Hematoxylin Solution for 1-2min, then wash with distilled water.(see note 5)
6. (optional) Differentiate by Acidic Differentiation Solution for 2-5s, then wash with distilled water.
7. (optional) Blue with Scott Bluing Solution for 3 min and wash with distilled water for 3min.
8. Dehydrate by series of ethanol, transparent by xylene and seal with resinene.

Result

Glycogen, neutral mucus and various glycoproteins	Purplish Red
Acidic mucus (sulfated and carboxylated)	Blue
Acidic polysaccharide, hyaluronic acid	Blue
Note: the cells or tissues that containing neutral mucus and acidic mucus can be dyed different shades of blue purple to purple.	

Note

1. Section dewaxing should be as clean as possible, otherwise it will affect the dyeing effect.
2. Reagent (B) has oxidation-reduction properties, while Reagent (C) has a pungent odor. It is recommended to take an appropriate amount, seal it in a timely manner, and store it at 2-8 °C.
3. The oxidation time of Oxidant should not be too long, and the best temperature is 18-22 °C.
4. After Schiff staining, the reagent should turn slightly red, which is a normal phenomenon. The tissue background should be colorless or light red.
5. The hematoxylin staining step is nuclear contrast staining, which can be omitted if the morphology is clear or HE control staining has been performed.
6. For your safety and health, please wear experimental clothes and disposable gloves.

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Alcian Blue Periodic Acid Schiff (AB-PAS) Stain Kit