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**Cat:** G1580

**Size:** 100T

**Storage:** -20°C, avoid light, valid for 1 year.

### Kit Components

Reagent	100T	Storage
Reagent(A): $\beta$ -Gal Fixative	100mL	2-8°C, avoid light
Reagent(B): X-Gal Solution	5mL	-20°C, avoid light
Reagent(C): $\beta$ -Gal Stain Solution A	1mL	2-8°C, avoid light
Reagent(D): $\beta$ -Gal Stain Solution B	1mL	2-8°C, avoid light
Reagent(E): $\beta$ -Gal Stain Solution C	100mL	2-8°C
Mix Reagent B,C,D,E, in proportion of 5:1:1:93 as Dyeing Working Solution. It is ready to use.		

### Introduction

The Senescence  $\beta$ -Galactosidase Staining Kit is a staining kit for senile cells or tissues based on the up-regulation of SA- $\beta$ -Gal(senescence-associated  $\beta$ -galactosidase) activity during aging. Aging of cells or tissues can be observed under optical microscopes. Cell aging is also considered as a way for organisms to inhibit tumors and a potential cause of organism aging. The kit can be used for senescence detection of cultured cells and tissue sections.

Using X-Gal as substrate, the Senescence  $\beta$ -Galactosidase Staining Kit produces dark blue products catalyzed by aging-specific  $\beta$ -galactosidase. Cells or tissues expressing  $\beta$ -galactosidase that turn blue are easily observed under optical microscopy. This kit only stains senescent cells. For tissue sections or blocks, the number of samples that can be detected depends on the size of the sample. For ordinary sections or 6-well plate, at least 100 samples can be detected.

### Protocol(for reference only)

*Note: It is recommended to use polypropylene (PP) container to prepare dyeing working solution, but polystyrene (PS) container cannot be used to prepare dyeing working solution.*

*The incubation at 37 °C during dyeing could not be carried out in carbon dioxide incubator.*

#### Adherent Cells Stain

- For cells cultured in 6-well plates, remove the cell culture medium and wash once with PBS or HBSS. Add 1 ml  $\beta$ -Gal Fixative to each pore and fix at room temperature for 15 min. For other types of culture plates, the amount of fixative solution and subsequent solution is operated according to this ratio.
- Remove  $\beta$ -Gal Fixative and wash cells 3 times with 1×PBS for 3 min each time. (*priority step see note 3*)
- Remove washing solution and add 1 ml Dyeing Working Solution to each hole.
- When incubate overnight for 12-24h at 37°C, the 6-well plate can be sealed with parafilm to prevent evaporation.
- View under the optical microscope. If you can't view and count in time, you can remove the Dyeing Working Solution, add 2ml PBS, and store it at 2-8 °C for several days; moreover, it can be stored for a long time after sealing with aqueous gelatin sealing agent.

#### Suspension Cell Stain

- Collect the cells into 1.5ml eppendorf tube by centrifugation, wash once with PBS or HBSS, and fix at room temperature for 15 min with 1 ml  $\beta$ -Gal Fixative. Slowly shake on the shaking table to avoid forming cell clumps.
- Centrifuge and remove  $\beta$ -Gal Fixative, wash cells with 1×PBS 3 times for 3 min each time. (*priority step see note 3*)
- Centrifuge and remove the washing solution then add 0.5-1 ml Dyeing Working Solution to each tube. Incubate overnight for 12-24h at 37°C.
- Drip some stained cells onto section and view under optical microscope. If not view in time, centrifuge and remove the dyeing working solution then adding 1 ml PBS at 2-8°C can keep for several days. If centrifuged, cells can also be taken for smears, and sealed with water based sealing agent, it can be stored for a long time at 2-8°C.

### Tissue Sections Stain

1. For enzymatically protected frozen slices, follow the steps below. (See Note 4)
2. Add a suitable volume of  $\beta$ -Gal Fixative to cover the tissue adequately and fix it for more than 15 min. (*priority step see note 3*)
3. Wash the tissues 3 times with PBS for more than 5 min each.
4. Remove washing solution.
5. When incubate overnight for 12-24h at 37°C, if drip dyeing, it can be placed in the wet-box to prevent evaporation; if soak dyeing, soak the section in dye container for closed condition.
6. View under optical microscope. If not view in time, can seal with glycerol gelatin for storage.

### Note

1. The  $\beta$ -Gal Fixative has certain corrosiveness and toxicity. Pay attention to the protection during operation.
2. The  $\beta$ -galactosidase staining reaction depends on specific pH conditions and cannot be performed in a carbon dioxide incubator. The high concentration of carbon dioxide in the carbon dioxide incubator used for cell culture can affect the pH value of the staining solution, leading to staining failure.
3. The priority step is using reagent E:  $\beta$ -Gal Stain Solution C to wash and incubate instead of 1×PBS before staining.
4. The stability of galactosidase is poor, and it is easy to inactivate when using paraformaldehyde to fix or prepare paraffin sections. It is recommended to prepare quick-frozen sections, or use enzyme protection special fixing solution: G2195- $\beta$ -Galactosidase Colorimetric Solution for fixing and freezing sections. Tissue immersion in galactosidase fixative stabilizes enzyme activity for at least one week.
5. PP containers or glass containers should be used in the preparation of dyeing working solution instead of PS containers. But dyeing can be carried out in PS containers, such as 6-well plate.
6. If the climbing film or smear needs to be stored for a long time, it can be washed twice with distilled water after dyeing, and the neutral gum seal can be made after the excess liquid is poured out and dried.
7. For your safety and health, please wear experimental clothes and disposable gloves.