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- [1] Xianghua Liu,Ning Sun,Nan Mo,et al. Quercetin inhibits kidney fibrosis and the epithelial to mesenchymal transition of the renal tubular system involving suppression of the Sonic Hedgehog signaling pathway. Food & Function. May 2019. (IF 3.289)
- [2] Houqi Ning,Yingqiu Li,Qiwen Tian,et al. The apoptosis of Staphylococcus aureus induced by glycinin basic peptide through ROS oxidative stress response. Food Science and Technology. January 2019;62-68. (IF 3.129)
- [3] Yong Zhang,Qian Sun,Ziwen Li,et al. Fermented soybean powder containing Bacillus subtilis SJLH001 protects against obesity in mice by improving transport function and inhibiting angiogenesis. Journal of Functional Foods. August 2019; 59:60-70. (IF 3.197)

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Cat: P1127

Size: 100mL

Storage: 2-8°C, avoid light, valid for 6 months.

Production

The purpose of immobilization is to preserve the original morphological structure of cells and tissues. Immobilizers can prevent endogenous lysosomal enzymes from autolysis of their own tissues and cells, and inhibit the growth of bacteria and molds. Fixed liquids are mainly divided into aldehyde based fixed liquids, mercury based fixed liquids, alcohol based fixed liquids, oxidant based fixed liquids, picrate based fixed liquids, etc. The more commonly used are aldehyde based and alcohol based fixed liquids.

Glutaraldehyde,4%(EM Grade) can fix more structures and chemical components, so the fixed components are not easily lost during the sample preparation process and will remain in place, making it suitable for electron microscopy observation and sample fixation. When observed under a light microscope, the coloring of each structure is generally darker and difficult to fade, resulting in lower recognition compared to conventional fixative solutions. Moreover, due to the stable fixation effect, dense tissue after fixation, and higher antigen masking ability than conventional fixation solution, appropriate adjustment of staining steps is required for routine staining.

Protocols (only for reference)

1. Operate according to the specific requirements of the experiment.
2. Take fresh specimens, wash them slightly with physiological saline, and then add Glutaraldehyde,4%(EM Grade) at 2-8 °C for 1-24 hours. For larger specimens, the fixation time should be appropriately extended.
3. Send for inspection or store at 2-8 °C.

Note

1. Gluta fixative has certain corrosiveness. Please handle it carefully in a well ventilated environment to avoid inhalation and skin contact.
2. The thickness of tissue sampling varies, and the fixation time also varies. It is recommended to take 2-4mm samples and fix them for at least 4 hours. Generally, it is not recommended to use a single sample exceeding 6mm. Proper material selection for the organization is beneficial for the penetration of the fixative.
3. The capacity of the fixed solution should be sufficient, and the volume ratio of the fixed solution to the tissue block should generally be greater than 10:1. If the volume is not large enough, the fixing solution can be replaced 1-3 times during the fixation period.
4. The temperature has a significant impact on fixation, and increasing the temperature can accelerate the fixation effect. However, it is best not to increase the temperature for this fixation solution.
5. For your safety and health, please wear laboratory clothes and disposable gloves for operation.