

GE Healthcare

Amersham High Molecular Weight Calibration Kit for SDS Electrophoresis

A lyophilized mixture of five highly purified well-characterized proteins for use in molecular weight determination in the presence of sodium dodecyl sulphate (SDS)

Product Booklet

Code: 17-0615-01



Page finder

1. Legal	3
2. Handling	4
2.1. Safety warnings and precautions	4
2.2. Storage	5
2.3. Expiry	5
3. Components	6
4. Other materials required	7
5. Description	8
6. Protocol	9
6.1. Preparation of calibration kit	9
6.2. Gel loading	10
6.3. Electrophoresis	10
6.4. Detection	10
6.5. Molecular weight determination	10
7. Typical results	12
8. Additional information	13
8.1. Background and references	13
9. Related products	14

1. Legal

GE and GE monogram are trademarks of General Electric Company. Amersham, CleanGel, ExcelGel, Mighty Small, Hoefer, Multiphor, PhastGel, PhastSystem and PlusOne are trademarks of GE Healthcare companies.

Coomassie is a trademark of Imperial Chemical Industries Ltd

© 2006 General Electric Company – All rights reserved.

GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specification and features shown herein, or discontinue the product described at any time without notice or obligation.

Contact your GE Healthcare representative for the most current information and a copy of the terms and conditions.

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited.

Amersham Place, Little Chalfont,
Buckinghamshire, HP7 9NA UK

2. Handling

2.1. Safety warnings and precautions

Warning: For research use only. Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

Human blood products provided as components of this pack have been obtained from donors who were tested individually and were found to be negative for the presence of Human Immunodeficiency Virus antibody (HIV-Ab)* as well as for Hepatitis B surface Antigen (HBsAg) using approved methods (ELISA)

As no test method can offer complete assurance that Hepatitis B virus, Human Immunodeficiency Virus antibody (HIV-Ab)* or other infectious agents are absent, all human blood products should be considered potentially infectious. Handling, use,

storage and disposal should be in accordance with the procedures defined by an appropriate National biohazard safety guideline or regulation, where it exists (for example USA Centre for Disease Control/National Institutes of Health manual "Biosafety in microbiological and Biomedical Laboratories", 2nd Edition 1988).

All chemicals should be considered as being potentially hazardous. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing, such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes, wash immediately with water. See material safety

data sheet(s) and/or safety statements for specific advice.

* HIV is the abbreviation used for HTLV-III and LAV.

2.2. Storage

The kit should be stored at 2–8°C.

2.3. Expiry

For expiry details see outer packaging.

3. Components

Protein mixture

176 µg/vial, 10 vials, each containing the following proteins:

Myosin (1), rabbit muscle, 25 µg, molecular weight (M_r) 220 000

α_2 -Macroglobulin (2), bovine plasma, 100 µg,

M_r 170 000

β -Galactosidase (3), *E.coli*, 16 µg, M_r 116 000

Transferrin (4), human,

17 µg, M_r 76 000

Glutamic dehydrogenase (5), bovine liver, 18 µg,

M_r 53 000

The amount of each protein has been chosen to give bands of equal intensity when stained with Coomassie™ Brilliant Blue following Laemmli-type gel electrophoresis. Intensities may vary when using other staining methods.

4. Other materials required

- Electrophoresis reagents appropriate to the application being run
- Detection reagents appropriate to the application being run
- Gel electrophoresis equipment

5. Description

The High Molecular Weight SDS Calibration Kit for SDS electrophoresis is a lyophilized mixture of five highly purified well-characterized proteins for use in molecular weight estimation in the presence of the detergent sodium dodecyl sulphate (SDS). The molecular mass of the protein under investigation is determined by comparing its electrophoretic mobility with that of proteins contained in the kit.

Ten vials are supplied, each containing a lyophilized mixture of highly purified protein standards of molecular mass range (M_r) 53 000 to 220 000 when used in denaturing polyacrylamide electrophoresis.

6. Protocol

6.1. Preparation of calibration kit

Reconstituting the contents of one vial in 100 µl of water gives a protein solution in 10 mM Tris-HCl (pH 8.0), 1 mM EDTA, 33 mM KCl, 2.5% dithiothreitol (DTT) and 2.5% SDS. The solution also contains 12% mannitol (stabilizer and density enhancer) and xylene cyanol green (tracking dye). It is not necessary to heat the calibration kit components in order to denature them before use. For best reproducibility, discard any unused portion of the reconstituted protein solution. However, if necessary, the solution can be stored at -80°C for 3 months.

For Coomassie Brilliant Blue detection

Aliquots from the reconstituted solution can be applied directly to the gel of choice for Coomassie Blue staining (Figure 2). However, if further dilutions are desired, use a standard 1x sample buffer (0.0625 M Tris-HCl, 2% SDS, 10% v/v glycerol, 0.1 M DTT and 0.01% bromophenol blue, pH 6.8).

For PhastGel™, ExcelGel™ and CleanGel™ precast gels, reconstitute the contents of a vial in 100 µl of 10 mM Tris-HCl, 2% SDS, 0.1 M DTT, 0.01% bromophenol blue and 1mM EDTA, pH 8.0.

For silver stain detection

For silver staining (Figure 3), reconstitute the contents of a vial as described for Coomassie blue staining, then dilute aliquots by at least 50-fold in 1 x sample buffer.

6.2. Gel loading

Select the appropriate sample volume from the table:

Gel type	sample volume (µl)
Vertical mini	3–10
Vertical standard	3–10
Multiphor™ II flatbed	2–4
PhastSystem™	0.3–4

6.3. Electrophoresis

Perform electrophoresis according to the instructions supplied with the gel apparatus being used.

6.4. Detection

Stain the gel using the desired method.

6.5. Molecular weight determination

Measure the migration distance of the proteins in the Calibration Kit and of the protein(s) of interest. Measure the migration distance of the dye marker. Calculate the corresponding R_f values by dividing migration distance of the protein by migration distance of the dye marker.

Construct a calibration curve by graphing R_f vs. log molecular weight for the proteins in the Calibration Kit (Figure 1). Determine the molecular weight of the protein(s) of interest from the calibration curve.




Figure 1. Calibration curve constructed using results shown in figure 2.

7. Typical results




Figure 2. HMW-SDS standards stained with Coomassie Brilliant Blue. Aliquots (10 μ l per lane) of a 2 x dilution were loaded on a self-cast 7.5% T, 2.7% C gel. The gel was run at a constant current of 20 mA for 1 hour, 42 minutes on a Mighty Small™ electrophoresis unit. The gel was stained with PhastGel Blue R (17-0518-01).

Protein	M_r (Da)	R_f
Myosin	220 000	0.02
α_2 -Macroglobulin	170 000	0.09
β -Galactosidase	116 000	0.24
Transferrin	76 000	0.41
Glutamic dehydrogenase	53 000	0.66




Figure 3. HMW-SDS standards stained with silver stain. Aliquots (5 μ l per lane) of a dilution series were loaded on an ExcelGel SDS Homogeneous 12.5(80-1261-01), run at 600 V, 50 mA, 30 W for 80 minutes on a Multiphor II flatbed unit. The gel was stained with PlusOne™ Silver Staining Kit, Protein (17-1150-01). The dilution factor, with respect to reconstitution of a vial in 100 μ l, is indicated in each lane.

Protein	M_r (Da)	R_f
Myosin	220 000	0.08
α_2 -Macroglobulin	170 000	0.12
β -Galactosidase	116 000	0.19
Transferrin	76 000	0.28
Glutamic dehydrogenase	53 000	0.40

8. Additional information

8.1. Background and references

For further information regarding molecular weight determinations and denaturing electrophoresis, see Hoefer Protein Electrophoresis Applications Guide (80-6013-88)

1. Starr, R., Dept of Clinical Veterinary Science, University of Bristol, UK, unpublished data.
2. Nakamura, S. et al., Hoppe-Seyler's Z. Physiol. Chem. **356**, 677-692 (1975).
3. Fowler, A. V. and Zabin, I., Proc. Natl. Acad. Sci. USA. **74**, 1507-1510 (1977)
4. Roberts, R. C. et al., J. Biol. Chem. **241**, 4907-4913 (1966).
5. Eisenberg, H. and Tomkins, G. M., J. Mol. Biol. **31**, 37-49, (1968).
6. Nellis, L. P. et al., Biochim. Biophys. Acta **623**, 46, (1980).

9. Related products

PhastGel Blue R (40 Coomassie Blue R-350 tablets)	17-0518-01
PlusOne silver Staining Kit, protein	17-1150-01
Hoefer™ Automated Gel Stainer with 19 x 29 cm PTFE coated staining tray	80-6395-02
with 29 x 35 cm PTFE coated staining tray	80-6396-16
Hoefer Protein Electrophoresis Application Guide	80-6013-88

GE Healthcare offices:

GE Healthcare Bio-Sciences AB
Björkgatan 30 751 84

Uppsala
Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5 D-79111
Freiburg
Germany

GE Healthcare UK Limited
Amersham Place
Little Chalfont
Buckinghamshire
HP7 9NA
UK

GE Healthcare Bio-Sciences
Corp
800 Centennial Avenue
P.O. Box 1327
Piscataway
NJ 08855-1327
USA

GE Healthcare Bio-Sciences KK
Sanken Bldg. 3-25-1
Hyakunincho Shinjuku-ku
Tokyo 169-0073
Japan

**GE Healthcare
regional office
contact numbers:**

Asia Pacific
Tel: +85 65 6 275 1830
Fax: +85 65 6 275 1829

Australasia
Tel: +61 2 8820 8299
Fax: +61 2 8820 8200

Austria
Tel: 01 /57606 1613
Fax: 01 /57606 1614

Belgium
Tel: 0800 73 890
Fax: 02 416 82 06

Canada
Tel: 1 800 463 5800
Fax: 1 800 567 1008

**Central, East, & South
East Europe**
Tel: +43 1 972720
Fax: +43 1 97272 2750

Denmark
Tel: 45 70 25 24 50
Fax: 45 16 24 24

Eire
Tel: 1 800 709992
Fax: 0044 1494 542010

Finland & Baltics
Tel: +358-(0)9-512 39 40
Fax: +358 (0)9 512 39 439

France
Tel: 01 6935 6700
Fax: 01 6941 9677

Germany
Tel: 0800 9080 711
Fax: 0800 9080 712

Italy
Tel: 02 26001 320
Fax: 02 26001 399

Japan
Tel: +81 3 5331 9336
Fax: +81 3 5331 9370

Korea
Tel: 82 2 6201 3700
Fax: 82 2 6201 3803

Latin America
Tel: +55 11 3933 7300
Fax: +55 11 3933 7304

Middle East & Africa
Tel: +30 210 9600 687
Fax: +30 210 9600 693

Netherlands
Tel: 0800 82 82 82 1
Fax: 0800 82 82 82 4

Norway
Tel: +47 815 65 777
Fax: 47 815 65 666

Portugal
Tel: 21 417 7035
Fax: 21 417 3184

**Russia & other C.I.S.
& N.I.S.**
Tel: +7 (495) 956 5177
Fax: +7 (495) 956 5176

South East Asia
Tel: 60 3 8024 2080
Fax: 60 3 8024 2090

Spain
Tel: 902 11 72 65
Fax: 935 94 49 65

Sweden
Tel: 018 612 1900
Fax: 018 612 1910

Switzerland
Tel: 0848 8028 10
Fax: 0848 8028 11

UK
Tel: 0800 515 313
Fax: 0800 616 927

USA
Tel: +1 800 526 3593
Fax: +1 877 295 8102

<http://www.gehealthcare.com/lifesciences>

GE Healthcare UK Limited

Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA
UK



imagination at work