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[1] Hongying Diao, Bin Liu, Yongfeng Shi, et al. MicroRNA-210 alleviates oxidative stress-associated cardiomyocyte apoptosis by regulating BNIP3. Bioscience. June 2017. (IF 5.876)

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Cat: DA1016

Size: 60mL/200mL

Storage: 2-8°C, avoid light, valid for 1 year.

Kit Components

Reagent	60mL	200mL	Storage
Reagent (A):Solution A	3mL	10mL	2-8°C, avoid light
Reagent (B):Solution B	3mL	10mL	2-8°C, avoid light
Reagent (C):Solution C	60mL	100mL×2	2-8°C
Before use, mix Reagent A, B and C as the ratio of 1:1:18 to prepare DAB Working Solution.			

Production

Enhanced DAB Chromogenic kit (20×) is a kind of reagent kit for immunohistochemical chromogenic, in situ hybridization chromogenic or Western, Southern, Northern, EMSA and other membrane chromogenic with the help of Horseradish peroxidase (HRP). DAB is a common substrate of Horseradish peroxidase. Under the catalysis of Horseradish peroxidase, DAB will produce brown precipitate. The brown precipitate is insoluble in water and ethanol and can be used for normal dehydration and sealing.

This product adopts a special formula with high sensitivity, low background, good repeatability, stable storage, and convenient use. It is suitable for staining and color reaction in protein blotting, immunohistochemistry, immunocytochemistry, dot blotting, and biochip.

Protocols (only for reference)

1. For tissue sections or protein imprinted membranes, after incubation with antibodies labeled with Horseradish peroxidase (HRP) or other forms of probes, rinse 3-5 times with appropriate washing solution for 3-5 minutes each time. (See Note 1)
2. Mix Solution A, B, and C evenly in a ratio of 1:1:18 to form DAB Working Solution. It is recommended to use within 4 hours. (See Note 2)
3. Add an appropriate amount of DAB Working Solution to the tissue slice or membrane to ensure sufficient coverage of the sample. (See Note 3)
4. Incubate at room temperature in dark for 1-30 minutes. If the chromogenic time is too long, it can cause an increase in background. Therefore, the chromogenic process should be closely observed (usually 3-10 minutes is optimal), and the chromogenic reaction should be terminated by rinsing with running water when the background is shallow and reaches the appropriate chromogenic intensity.
5. For tissue sections or cell samples, after the chromogenic reaction is terminated, other dyes (commonly used hematoxylin staining solution) can be used for re staining. For the film, after the chromogenic reaction is terminated, it can be dried at room temperature and stored away from light.

Note

1. Rinsing solution can use 1×PBS, 1×PBST and other commonly used washing buffers in immunohistochemistry or membrane chromogenic.
2. Reagent (A): Solution A should be a light yellow to purplish red transparent liquid, and the newly prepared working solution should be colorless or brownish yellow. If the color is too dark or precipitates appear, do not use it. If the prepared working solution has been left for a long time, causing discoloration or sedimentation, please do not use it.
3. It is usually recommended to add 80-100ul DAB Working Solution to each tissue section, and 40-50ul DAB Working Solution to each square centimeter of imprinting film. It is that a 60mL Enhanced DAB Chromogenic kit (20×) can be used for staining approximately 600 tissue sections or staining imprinted membranes with an area of 1200 square centimeters.
4. DAB has potential mutagenic effect. For your safety and health, please wear White coat and disposable gloves.
5. If the chromogenic reaction is poor or cannot confirm the reagent quality due to prolonged non-use of reagents, can mix diluted secondary antibody incubation solution and DAB Working Solution in a ratio of 1:3, and left to stand for observation. The mixture should turn brownish yellow within 2 minutes and precipitate within 10 minutes.