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## Peroxidase Blocking Buffer, Enhanced

**Cat:** G2850

**Size:** 100mL

**Storage:** 2-8°C, avoid light, valid for 1 year.

### Introduction

Endogenous peroxidase is widely present in some cells and tissues, including red blood cells, kidneys, and liver tissues. It is prone to high background and even false positive results during HRP detection system coloration. Therefore, it is advisable to use appropriate peroxidase blocking solution for cell or tissue samples before staining to eliminate interference from endogenous peroxidase.

The effectiveness of the peroxidase strong blocking solution produced by Solebao is significantly superior to conventional hydrogen peroxide solution and inhibitor type blocking solution, which can fully reduce non-specific background staining and have a certain enhancement effect on specific staining. It can be used for blocking endogenous peroxidase in tissues or cells during immunohistochemistry (IHC), immunocytochemical staining (ICC), and in situ hybridization (ISH), and can also be used for step-by-step masking in tyramide fluorescence multi labeling.

### Protocol(for reference only)

#### For IHC and ICC Sample

1. For samples subjected to IHC and ICC staining, Wash the sample with P1031-1×PBST buffer solution or a liquid with cell membrane permeability function for 3-5 min first.
2. Absorb the washing solution thoroughly, add an appropriate amount of Peroxidase Blocking Buffer, Enhanced dropwise to completely cover the sample, and incubate at room temperature for 5-10 min or 37 °C for 2-5 min.
3. Remove the Peroxidase Blocking Buffer, Enhanced, and wash the sample 2-3 times with washing solution or PBS or other appropriate solutions, each time for 2-5 min. Then subsequent steps such as site closure and primary antibody incubation can be carried out.

#### For ISH Sample

1. Before conducting probe hybridization, add an appropriate amount of Peroxidase Blocking Buffer, Enhanced dropwise to cover the sample, and incubate at room temperature for 5-10 min or 37°C for 2-5 min.
2. Remove the blocking solution and wash 2-3 times with the washing solution of in situ hybridization, each time for 2-5 min. Subsequently, subsequent steps such as hybridization can be carried out.

### Note

1. This product has tested the immunostaining effects of various antibodies, which are equivalent to or significantly improved by conventional methods.
2. For some special situations where it is difficult to use a peroxidase detection system, it is possible to consider using an alkaline phosphatase detection system.
3. This product does not contain highly toxic methanol, but it is still irritating to the human body. Please pay attention to appropriate protection.
4. This product is limited to scientific research by professionals and shall not be used for clinical diagnosis or treatment. It shall not be used for food or medicine, and shall not be stored in ordinary residential areas.
5. For your safety and health, please wear laboratory clothes and disposable gloves for operation.