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# Tyramine Signal Amplification Kit, with HRP-Goat anti-Rabbit IgG<sup>V02</sup> and TYR-Cy5

**Cat:** G4867

**Size:** 100T/300T

**Storage:** -20°C, avoid light, valid for 1 year. Period after opening is 3 months.

## Kit Components

| Reagent                                                                                       |                                       | 100T  | 300T  | Storage            |
|-----------------------------------------------------------------------------------------------|---------------------------------------|-------|-------|--------------------|
| Reagent(A):<br>Tyr-label<br>Working Solution                                                  | Reagent(A1):TYR-Cy5 Solution 500×     | 20uL  | 60uL  | -20°C, avoid light |
|                                                                                               | Reagent(A2):Enhancement Solution 100× | 0.1mL | 0.3mL | 2-8°C, avoid light |
|                                                                                               | Reagent(A3):Diluent Solution          | 10mL  | 30mL  | 2-8°C              |
| Mix A1: A2: A3 in the ratio of 1:5:500 to prepare Tyr-label Working Solution, and use it now. |                                       |       |       |                    |
| Reagent(B):Peroxidase Enhanced Blocking Solution                                              |                                       | 10mL  | 30ml  | RT                 |
| Reagent(C):HRP-Goat anti-Rabbit IgG 100×                                                      |                                       | 0.1mL | 0.3mL | -20°C, avoid light |

## Introduction

The tyramide signal amplification (TSA) system can be used to detect low abundance targets in fluorescence immunocytochemistry (ICC), immunohistochemistry (IHC), and in situ hybridization (ISH) techniques, which can increase signal sensitivity by 100 times. The TSA fluorescence kit uses horseradish peroxidase (HRP) to directly catalyze the covalent deposition of fluorescent groups around the immobilized enzyme. The labeling process is rapid (less than 10 min), and the deposition labeling can be directly observed under a standard or confocal microscope. The TSA kit can be combined with traditional staining methods for multicolor imaging, or more tyramine reactions can be sequentially performed to label different targets on a sample.

The Solarbio Tyramide Signal Amplification (TSA) series products includes 8 kits: the G4860/G4861 is fluorescein labeled and detected at 494nm/ 517nm; The G4862/G4863 was fluorescently labeled as TRAMA and detected at 552nm/ 578nm; The G4864/G4865 was fluorescently labeled as Cyanine3 and detected at 550nm/ 570nm; The G4866/G4867 was fluorescently labeled as Cyanine5 and detected at 648nm/667nm.

## Self Provided Materials

1×PBS, 4% PFA, Translucent (Triton X-100, strong), 1×PBST, Sodium Citrate Antigen Repair Solution (1×) or EDTA Antigen Repair Solution, 5% BSA Blocking Solution, or Serum-free Blocking Solution. The above reagents can be self prepared or purchased as commercially available finished reagents.

## Protocols(for reference only)

### Sample Preparation

- (Optional) Prepare a negative control sample (sample without incubating primary antibody). (See Note 7)
- Cell sample preparation:
  - ① Wash cells with 1×PBS twice for 2 min each time.
  - ② After removing PBS, add an appropriate amount of pre-cooled 4% tissue cell fixative and fix at RT for 15 min. (See Note 1)
  - ③ Wash cells with 1×PBS twice for 2 min each time.
- Preparation of paraffin sections:
  - ① Place the paraffin slices in a 60 °C oven for 30 min.
  - ② Soak in xylene or Tissue Transparent Dewaxing Solution, Benzene Free twice at RT for 5 min each time to thoroughly dewaxing.
  - ③ At RT, the series of ethanol rehydration:100% ethanol I, 100% ethanol II, 95% ethanol, 90% ethanol, 80% ethanol, 70% ethanol, and distilled water for 5 min each. Soak in 1×PBS for 3 min.
- Preparation of frozen sections:
  - ① Remove the frozen slices from the refrigerator and place them in a 2-8 °C refrigerator for 10 min to balance. Wash sections with 1×PBS twice for 2 min each time.
  - ② After absorbing PBS, fix with pre-cooled 4% PFA and fix at RT for 15 min. (See Note 2)
  - ③ Wash the sample twice with 1×PBS for 2 min each time.
- Use a pap pen to outline the sample for subsequent localization and imaging.

### Repair and Closure

1. (Optional) Antigen repair: Heat the antigen repair solution adapted to the first antibody to boiling, place the prepared tablets in buffer, and intermittently boil for 10 min. After antigen repair, place the repair solution together at RT and slowly cool down. (See Note 3)
2. (Optional) Cell permeability: Add an appropriate amount of permeability agent (Triton X-100, strong) to the sample and incubate at RT for 10 min. (See Note 4)
3. (Optional) Peroxidase inactivation: Add an appropriate amount of reagent to each sample (B): Peroxidase Enhanced Blocking Solution at RT for 5-10 min, wash cells with 1×PBS twice for 2 min each time.
4. Site blocking: Cover the sample at RT and block it for 10min using serum or BSA or serum-free blocking solution.

### Immuno Labeling

1. Dilute the first antibody with sealing solution, dropwise coverage sample and incubate at 37 °C for 1 h or 2-8 °C overnight (8-12 h).
2. Wash with 1×PBS at RT 3 times for 5 min each time.
3. Dilute the Reagent(C):HRP-Goat anti-Rabbit IgG 100× with blocking solution, dropwise coverage sample and incubate at 37 °C for 40 min.
4. Wash with 1×PBS at RT 3 times for 5 min each time.
5. Prepare 100μL-300μL Tyr-label Working Solution for each sample. The working fluid is stored at RT in a dark place and can be stored overnight at most.
6. Cover the sample with dropwise incubation solution and incubate it in dark at RT for 10 min.
7. Wash with 1×PBS at RT 3 times for 5 min each time.
8. Microscopic imaging. For tissue samples on glass slides, please use water-based gelatin sealing agent or anti fluorescence attenuation sealing agent before microscopic imaging.

### Note

1. The fixed solution and dewaxing reagents involved in the operation steps have certain risks, and it is recommended to operate in a fume hood.
2. Calculate the reagent dosage for each slice based on the tissue area, usually 100uL/cm<sup>2</sup>, and add half of the reference medium volume to the culture plate.
3. Cross linked fixatives such as formaldehyde and glutaraldehyde can stabilize tissue structure but also mask antigen binding sites to a certain extent. Therefore, it is recommended to perform antigen repair on conventional fixed samples. If ethanol acetone is used for fixation, it can be omitted.
4. A complete membrane structure can hinder the entry and exit of antibodies and fluorescent dyes, so it is recommended to conduct transparent treatment for frozen sections and cell samples.
5. Please operate different fluorescence signal amplification kits according to the recommended wavelength.
6. Compared with fluorescent secondary antibodies, the tyramide signal amplification kit can display higher sensitivity and signal intensity. Therefore, during the experiment, the concentration of the primary antibody can be appropriately reduced to reduce the background fluorescence caused by non-specific binding. We suggest setting a concentration gradient of the primary antibody to find the optimal dilution ratio.
7. If background fluorescence needs to be considered, it is recommended to set a negative control that has not been incubated with the first antibody. Ensure that the negative control is not cross contaminated by reagents in the positive sample during incubation and washing processes. For tissue samples, we also recommend imaging the unstained control (without adding antibodies or tyramides) to determine the effect of tissue spontaneous fluorescence on the background.
8. A high concentration of tyramine substrate may lead to strong signal or high background, and the best imaging effect can be obtained by searching from 1:5:200 to 1:5:1000.
9. Multiple tyramide amplification kits can be used in sequence to label different targets on the same sample by HRP quenching or antibody stripping after each tyramide reaction.
10. For your safety and health, please wear laboratory clothes and disposable gloves for operation.
11. This product is for scientific research purposes only. Do not use for medical, clinical, food, or cosmetic purposes.