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Four Color Tyramine Signal Amplification Kit, with HRP-Goat^{V02} anti-Rabbit IgG

Cat: G4881

Size: 100T/300T

Storage: -20°C, avoid light, valid for 1 year. Period after opening is 3 months.

Kit Components

Reagent		100T	300T	Storage
Reagent(A):TYR-FITC Solution 500×		20uL	60uL	-20°C, avoid light
Reagent(B):TYR-CY3 Solution 500×		20uL	60uL	-20°C, avoid light
Reagent(C):TYR-CY5 Solution 500×		20uL	60uL	-20°C, avoid light
Reagent(D): Tyr-label Diluent	Reagent(D1):Enhancement Solution 100×	0.4mL	1mL	2-8°C, avoid light
	Reagent(D2):Diluent Solution	40mL	100mL	2-8°C
Mix D1,D2 by 1: 100 to prepare Tyr-label Diluent.Then add A or B or C in Tyr-label Diluent by 1:500 to form Tyr-label Working Solution.				
Reagent(E):Peroxidase Enhanced Blocking Solution		40mL	100mL	RT
Reagent(F):HRP-Goat anti-Rabbit IgG 100×		0.4mL	1mL	-20°C, avoid light
Reagent(G):Special Antifade Mountant		20mL	60mL	-20°C, avoid light

Introduction

The tyramide signal amplification (TSA) system can be used to detect low abundance targets in fluorescence immunocytochemistry (ICC), immunohistochemistry (IHC), and in situ hybridization (ISH) techniques, which can increase signal sensitivity by 100 times. The TSA fluorescence kit uses horseradish peroxidase (HRP) to directly catalyze the covalent deposition of fluorescent groups around the immobilized enzyme. The labeling process is rapid (less than 10 min), and the deposition labeling can be directly observed under a standard or confocal microscope. The TSA kit can be combined with traditional staining methods for multicolor imaging, or more tyramine reactions can be sequentially performed to label different targets on a sample.

The Solarbio Polychromatic Tyramide Signal Amplification (TSA) series products includes 2 four color multi label fluorescence enhancement kits, labeled with fluorescein, Cyanine3, Cyanine5, and DAPI, which can simultaneously label three proteins and nuclei.

Self Provided Materials

1×PBS, 4% PFA, Translucent (Triton X-100, strong), 1×PBST, Sodium Citrate Antigen Repair Solution (1×) or EDTA Antigen Repair Solution, 5% BSA Blocking Solution, or Serum-free Blocking Solution. The above reagents can be self prepared or purchased as commercially available finished reagents.

Protocols(for reference only)

Sample Preparation

- (Optional) Prepare a negative control sample (sample without incubating primary antibody). (See Note 7)
- Cell sample preparation:
 - ① Wash cells with 1×PBS twice for 2 min each time.
 - ② After removing PBS, add an appropriate amount of pre-cooled 4% tissue cell fixative and fix at RT for 15 min.Wash cells with 1×PBS twice for 2 min each time. (See Note 1)
- Preparation of paraffin sections:
 - ① Place the paraffin slices in a 60 °C oven for 30 min.
 - ② Dewax by xylene or Tissue Dewaxing Solution, Benzene Free twice at RT for 5 min each time.
 - ③ At RT, the series of ethanol rehydration:100% ethanol I, 100% ethanol II, 95% ethanol, 90% ethanol, 80% ethanol, 70% ethanol, and distilled water for 5 min each.Soak in 1×PBS for 3 min.
- Preparation of frozen sections:
 - ① Remove the frozen slices from the refrigerator and place them in a 2-8 °C refrigerator for 10 min to balance. Wash sections with 1×PBS twice for 2 min each time.
 - ② After absorbing PBS, fix with pre-cooled 4% PFA and fix at RT for 15 min. (See Note 2)
 - ③ Wash the sample twice with 1×PBS for 2 min each time.

5. Use a pap pen to outline the sample for subsequent localization and imaging.

Repair and Closure

1. (Optional) Antigen repair: Heat the antigen repair solution adapted to the first antibody to boiling, place the prepared tablets in buffer, and intermittently boil for 10 min. After antigen repair, place the repair solution together at RT and slowly cool down. (See Note 3)
2. (Optional) Cell permeability: Add an appropriate amount of permeability agent (Triton X-100, strong) to the sample and incubate at RT for 10 min. (See Note 4)
3. (Optional) Peroxidase inactivation: Add an appropriate amount of reagent to each sample (E): Peroxidase Enhanced Blocking Solution at RT for 5-10 min, wash cells with 1×PBS twice for 2 min each time.
4. Site blocking: Cover the sample at RT and block it for 10min using serum or BSA or serum-free blocking solution.

Immune Labeling

1. Dilute the first antibody with sealing solution, dropwise coverage sample and incubate at 37 °C for 1 h or 2-8 °C overnight (8-12 h). Wash with 1×PBS at RT 3 times for 5 min each time.
2. Dilute the Reagent(F):HRP-Goat anti-Rabbit IgG 100× with blocking solution, dropwise coverage sample and incubate at 37 °C for 40 min. Wash with 1×PBS at RT 3 times for 5 min each time.
3. Prepare 100μL-300μL Tyr-label Working Solution for each sample. The working fluid is stored at RT in a dark place and can be stored overnight at most.
4. Cover the sample with dropwise incubation solution and incubate it in dark at RT for 10 min.
5. Wash with 1×PBS at RT 3 times for 5 min each time.
6. Microscopic imaging. Sealing with Reagent(G):Special Antifade Mountant before microscopic imaging.

Continuous Marking

1. Repeat the repair, Closure, and immune labeling processes, using Tyr-label Working Solution prepared with different primary antibody proteins of the same genus and different TYR fluorescence labeled reserve solutions for single slice multiple labeling.
2. Use Reagent(G):Special Antifade Mountant to seal the film and then perform microscopic imaging.

Note

1. The fixed solution and dewaxing reagents involved in the operation steps have certain risks, and it is recommended to operate in a fume hood.
2. Calculate the reagent dosage for each slice based on the tissue area, usually 100uL/cm², and add half of the reference medium volume to the culture plate.
3. Cross linked fixatives such as formaldehyde and glutaraldehyde can stabilize tissue structure but also mask antigen binding sites to a certain extent. Therefore, it is recommended to perform antigen repair on conventional fixed samples. If ethanol acetone is used for fixation, it can be omitted.
4. A complete membrane structure can hinder the entry and exit of antibodies and fluorescent dyes, so it is recommended to conduct transparent treatment for frozen sections and cell samples.
5. Please operate different fluorescence signal amplification kits according to the recommended wavelength.
6. Compared with fluorescent secondary antibodies, the tyramide signal amplification kit can display higher sensitivity and signal intensity. Therefore, during the experiment, the concentration of the primary antibody can be appropriately reduced to reduce the background fluorescence caused by non-specific binding. We suggest setting a concentration gradient of the primary antibody to find the optimal dilution ratio.
7. If background fluorescence needs to be considered, it is recommended to set a negative control. Ensure that the negative control is not cross contaminated by reagents in the positive sample during incubation and washing processes. For tissue samples, we also recommend imaging the unstained control (without adding antibodies or tyramides) to determine the effect of tissue spontaneous fluorescence on the background.
8. A high concentration of tyramine substrate may lead to strong signal or high background, and the best imaging effect can be obtained by searching from 1:5:200 to 1:5:1000.
9. Multiple tyramide amplification kits can be used in sequence to label different targets on the same sample by HRP quenching or antibody stripping after each tyramide reaction.
10. For your safety and health, please wear laboratory clothes and disposable gloves for operation.
11. This product is for scientific research only. Do not use for medical, clinical, food, or cosmetic purposes.