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## Malachite Green Aqueous Solution, 5%

**Cat:** G1630 **Size:** 100mL

**Storage:** RT, avoid light, valid for 6 months.

### Introduction

Malachite green is a synthetic dye with CAS No. 569-64-2 and molecular weight of 364.92.

Malachite Green Aqueous Solution, 5% is mainly composed of malachite green and deionized water, which can be used for Brucella staining or other purposes.

## **Protocol** (for reference only)

- 1. Malachite Green Solution is generally used for special purposes. If it is used for general dyeing, it should be diluted to 0.5-1% to form Malachite Green Working Solution.
- 2. Smear: number the slide, prepare bacterial smears, add 2-3 drops of Malachite Green Working Solution on the sample of the slide, and dye for 2-3mins.
- 3. Cover the clean cover glass, and press it gently to make a tablet.
- 4. View the morphology of bacteria by microscope under low, high or oil microscope.

#### Note

- 1. Before smearing, make circle mark on the back in advance to judge the position of subsequent test.
- 2. The culture time of bacteria will affect the staining. If the culture time of positive bacteria is too long or it has died or the bacteria is dissolved, it often shows negative reaction.