

硫酸链霉素说明书

货号: S8290

规格: 5g/25g

保存: 2-8°C

效价: 650-850mcg/mg

产品说明:

抗菌素，主要抗革兰氏阴性菌，常与青霉素一起用于细胞培养，防止微生物污染。使用浓度一般为 100µg/mL。易溶于水，10ug/ml 即可达到抑菌效果。

相关文献:

- [1] Lulu Xie,Hanhan Dai,Minjing Li,et al. MARCH1 encourages tumour progression of hepatocellular carcinoma via regulation of PI3K-AKT- β -catenin pathways. Journal of Cellular and Molecular Medicine. May 2019; 23 (5):3386-3401. (IF 4.658)
- [2] Dandan Sheng,Jinxu Li,Chengchong Ai,et al. Electrospun PCL/Gel aligned scaffolds enhance the biomechanical strength in tendon repair. Journal of Materials Chemistry B. June 2019. (IF 5.047)
- [3] Xinnan Liu,Wentao Zhang,Lunjie Huang,et al. Fluorometric determination of dopamine by using molybdenum disulfide quantum dots. Microchimica Acta. April 2018. (IF 5.705)
- [4] Huateng Zhou,Lixiang Wang,Cui Zhang,et al. Feasibility of repairing full-thickness skin defects by iPSC-derived epithelial stem cells seeded on a human acellular amniotic membrane. Stem Cell Research & Therapy. May 2019. (IF 4.627)

注: 更多使用本产品的文献请参考索莱宝官网。