

PG 固定液 (2%/2.5%)

货号: P1116

规格: 100mL /500mL

保存: 2-8℃避光保存, 有效期 1 年。

产品简介:

组织和细胞常用的固定液, 尤其适用在电镜中保护超微结构, 不适用于固定核酸、多糖物质。主要由 2%多聚甲醛、2.5%戊二醛、磷酸盐组成,pH 值为 7.4, 有刺鼻气味。

使用说明 (仅供参考):

离心收集细胞, PBS 漂洗两遍, 2%多聚甲醛-2.5%戊二醛固定液固定 30min-1h, PBS 漂洗两遍, 用于后续实验。

相关试剂:

P2210	Bouin 氏固定液
P2220	carnoy 固定液
C1010	0.01M 柠檬酸钠缓冲液 pH6.0
P1110	4%多聚甲醛
P8400	PBS 干粉 (0.01M, PH7.2-7.4)
SA0011	SABC(小鼠 IgG)-POD kit
C1033	EDTA 抗原修复液

注意事项:

1. 产品信息仅供参考, 如有疑问请致电 400-968-6088 咨询。
2. 本产品仅供科研使用。请勿用于医药、临床诊断或治疗, 食品及化妆品等用途。请勿存放于普通住宅区。
3. 为了您的安全和健康, 请穿好实验服并佩戴一次性手套和口罩操作。

相关文献:

- [1] Xianghua Liu,Ning Sun,Nan Mo,et al. Quercetin inhibits kidney fibrosis and the epithelial to mesenchymal transition of the renal tubular system involving suppression of the Sonic Hedgehog signaling pathway. Food & Function. May 2019. (IF 3.289)
- [2] Saqi Er,Shuyin Bao,Shuying Han,et al. Effects of Agriophyllum squarrosum extracts on glucose metabolism in KKAY mice and the associated underlying mechanisms. Journal of Ethnopharmacology. September 2019. (IF 3.414)
- [3] Yong Zhang,Qian Sun,Ziwen Li,et al. Fermented soybean powder containing Bacillus subtilis SJLH001 protects against obesity in mice by improving transport function and inhibiting angiogenesis. Journal of Functional Foods. August 2019; 59:60-70. (IF 3.197)
- [4] Xia Huang,Feilong Wang,Chen Zhao,et al. Dentinogenesis and Tooth-Alveolar Bone Complex Defects in BM P9/GDF2 Knockout Mice. Stem Cells and Development. May 2019;28(10). (IF 3.147)

注: 更多使用本产品的文献请参考索莱宝官网。

