



羧基化四氧化三铁纳米颗粒

货号: NM000610

规格: 10mg/20mg

保存: 4°C避光密封保存, 保质期为 6 个月

Sealed, avoid light, and keep at 4°C. Expiry date: 6 months.


产品说明:

名称 (Product name)	羧基化四氧化三铁纳米颗粒 (DMSA Coated Fe ₃ O ₄ Nanoparticles)
颗粒直径 (Size)	~7, 20, 30, 50nm (TEM)
颜色 (Color)	棕色 (Brown)
浓度 (Concentration)	~4 mg/ml (Fe)

备注: 颗粒直径为单次测量结果, 不同批次间允许少量波动

Note: The particle size is a single measurement result, and a small amount of deviation is allowed.

表征测试图:



羧基化四氧化三铁纳米颗粒(7nm)的 TEM 图
TEM image of DMSA Coated Fe₃O₄ Nanoparticles (7nm)

应用:

用于磁共振成像对比、生物分子偶联、固定, 纳米探针构建等。

Magnetic resonance imaging contrast, biomolecule coupling, immobilization, nano probe construction, etc.

注意事项:

1. 本产品仅供科研使用。请勿用于医药、临床诊断或治疗, 食品及化妆品等用途。请勿存放于普通住宅区。
2. 为了您的安全和健康, 请穿好实验服并佩戴一次性手套和口罩操作。