

A498 Whole Cell RIPA Extract

Cat: SCC-110511-L10

Summary:

【Product name】: A498 Whole Cell RIPA Extract

【Quantity size】: 100ug

【Concentration】: 2.5mg/ml

【Isotype】: SDS-PAGE Protein Sample; Cell Lysate

【Tested applications】: WB; SDS-PAGE

【Storage】: Store at -20°C. Avoid freeze / thaw cycles.

Background:

A498 cells were lysed using radio immunoprecipitation assay (RIPA) buffer supplemented with protease and phosphatase inhibitors, centrifuged and the supernatant collected. 100 μ g of cellular protein was mixed with 4 times volume of 5X SDS-PAGE loading buffer (with DTT).