

U251 Whole Cell RIPA Extract

Cat: SCC-111712-L10

Summary:

【Product name】 : U251 Whole Cell RIPA Extract

【Quantity size】 :100ug

【Concentration】 : 2.5mg/ml

【Isotype】 : SDS-PAGE Protein Sample;Cell Lysate

【Tested applications】 : WB;SDS-PAGE

【Storage】 : Store at -20°C. Avoid freeze / thaw cycles.

Background:

U251 cells were lysed using radio immunoprecipitation assay (RIPA) buffer supplemented with protease and phosphatase inhibitors, centrifuged and the supernatant collected.100 μ g of cellular protein was mixed with 4 times volume of 5X SDS-PAGE loading buffer(with DTT).