

Jurkat Whole Cell RIPA Extract

Cat: SCC-121815-L10

Summary:

【Product name】: Jurkat Whole Cell RIPA Extract

【Quantity size】:100ug

【Concentration】: 2.5mg/ml

【Isotype】: SDS-PAGE Protein Sample;Cell Lysate

【Tested applications】: WB;SDS-PAGE

【Storage】: Store at -20°C. Avoid freeze / thaw cycles.

Background:

Jurkat cells were lysed using radio immunoprecipitation assay (RIPA) buffer supplemented with protease and phosphatase inhibitors, centrifuged and the supernatant collected.100 μ g of cellular protein was mixed with 4 times volume of 5X SDS-PAGE loading buffer(with DTT).