

## Cy5.5 Conjugated Antibody

### Summary:

**Quantity size:** 0.1mL/0.5mL

**Concentration:** 2mg/mL Buffer=0.01M PBS (pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

**Storage:** Shipped at 4°C, Store at -20°C, at least one year (Avoid repeated freeze/thaw cycles).

**Purification method:** Protein A/G

### Background:

Cyanine dyes are fluorescent dyes with excellent properties. Cy5.5 conjugated antibody is a kind of immunochemistry reagents widely used in immunology, molecular biology and various branches of clinical medicine. The various Cy5.5-labeled antibodies developed by Solarbio are all self-prepared immune serums. The immune serums are purified by affinity chromatography to obtain high-purity IgG.

### Application:

IF=1:50-500

Excitation spectrum: 675nm

Emission spectrum: 691nm

Not yet tested in other applications.

Dilution is subject to the product label. Optimal working dilutions must be determined by the end user.

### Related products:

SPA131 Goat anti-mouse IgG (purified)

SA134 Goat anti-rabbit IgG (immune serum)


SF134 Goat anti-rabbit IgG-FITC

SE237 Rabbit anti-pig IgG-HRP

K1031G-Cy3 Goat anti-mouse IgG/ Cy3

K1034G-Cy5 Goat anti-rabbit IgG/ Cy5

### Product Image:



Immunofluorescence analysis of Goat anti-Mouse IgG/Cy5.5 (K1031G-Cy5.5) was performed using SHSY5Y cells stained with Anti-KRT19 Monoclonal Antibody (K200066M). The cells were fixed with 4% paraformaldehyde for 10 minutes, permeabilized with 0.1% Triton X-100 for 10 minutes, blocked with 1% BSA for 1 hour and labeled with Anti-KRT19 Monoclonal Antibody for 2 hours at room temperature. Goat anti-Mouse IgG/Cy5.5 (K1031G-Cy5.5) was used at 1: 400 in phosphate buffered saline containing 0.2% BSA for 45 minutes at room temperature, for detection of KRT19 in the cytoplasm (red). Nuclei (blue) were stained with DAPI (C0065).

**Important Note:** This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.