

Rabbit anti-Goat IgG (Immune serum)

Cat: SA238

Summary:

【Product name】 :Rabbit anti-Goat IgG (Immune serum)

【Type】 : Polyclonal Antibody

【Source】 : Rabbit

【Species reactivity】 : Goat

【Immunogen】 :Goat IgG full length protein

【Calculated】 : MW:150kDa

【Tested applications】 : ELISA WB Double agar diffusion assay

【Recommended dilution】 :WB 1:1000-20000.;ELISA:1:5000-400000.;Double agar diffusion assay 1:4-32.

【Product form】 :Liquid

【Storage】 : Store at -20℃ at least one year (Avoid repeated freeze/thaw cycles).

Background:

Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10¹¹ variants. This product is reddish transparent liquid, containing 0.03% Proclin300.