



DMEM 培养基

货号: 31600/11995/11965/12100/12110

规格: 500mL

保存: 2-8°C, 避光保存, 有效期 12 个月。

产品内容:

Cat. NO.	31600	11995	11965	12100	12110
Components	Concentration (mg/L)	Concentration (mg/L)	Concentration (mg/L)	Concentration (mg/L)	Concentration (mg/L)
Glycine	30.0	30.0	30.0	30.0	30.0
L-Arginine hydrochloride	84.0	84.0	84.0	84.0	84.0
L-Cystine•2HCl	63.0	63.0	63.0	63.0	63.0
L-Glutamine	584.0	584.0	584.0	584.0	584.0
L-Histidine hydrochloride•H ₂ O	42.0	42.0	42.0	42.0	42.0
L-Isoleucine	105.0	105.0	105.0	105.0	105.0
L-Leucine	105.0	105.0	105.0	105.0	105.0
L-Lysine hydrochloride	146.0	146.0	146.0	146.0	146.0
L-Methionine	30.0	30.0	30.0	30.0	30.0
L-Phenylalanine	66.0	66.0	66.0	66.0	66.0
L-Serine	42.0	42.0	42.0	42.0	42.0
L-Threonine	95.0	95.0	95.0	95.0	95.0
L-Tryptophan	16.0	16.0	16.0	16.0	16.0
L-Tyrosine disodium salt dihydrate	104.0	104.0	104.0	104.0	104.0
L-Valine	94.0	94.0	94.0	94.0	94.0
Choline chloride	4.0	4.0	4.0	4.0	4.0
D-Calcium pantothenate	4.0	4.0	4.0	4.0	4.0
Folic Acid	4.0	4.0	4.0	4.0	4.0
Niacinamide	4.0	4.0	4.0	4.0	4.0
Pyridoxine hydrochloride	4.0	4.0	4.0	4.0	4.0
Riboflavin	0.4	0.4	0.4	0.4	0.4
Thiamine hydrochloride	4.0	4.0	4.0	4.0	4.0
i-Inositol	7.2	7.2	7.2	7.2	7.2
Calcium Chloride (CaCl ₂) (anhyd.)	200.0	200.0	200.0	200.0	200.0
Ferric Nitrate (Fe(NO ₃) ₃ •9H ₂ O)	0.1	0.1	0.1	0.1	0.1
MgSO ₄ (anhyd.)	97.67	97.67	97.67	97.67	97.67
Potassium Chloride (KCl)	400.0	400.0	400.0	400.0	400.0
Sodium Bicarbonate (NaHCO ₃)	3700	3700.0	3700.0	3700	3700
Sodium Chloride (NaCl)	6400.0	6400.0	6400.0	6400.0	6400.0
NaH ₂ PO ₄ •H ₂ O	125.0	125.0	125.0	125.0	125.0
D-Glucose (Dextrose)	1000.0	4500.0	4500.0	4500.0	4500.0
Phenol Red	15.0	15.0	15.0	15.0	15.0
Sodium Pyruvate	110.0	110.0	-	-	-
HEPES	3600	-	-	3600	3600
Penicillin G Sodium salt	66.6	-	-	66.6	-
Streptomycin Sulfate	100	-	-	100	-

注: 更多使用本产品的文献请参考索莱宝官网。