



RPMI Medium 1640

货号: 10491、11875、31800、90023

规格: 500mL

保存: 2-8°C, 避光保存, 有效期 12 个月。

产品内容:

Cat. NO.	10491	11875	31800	90023
Components	Concentration (mg/L)	Concentration (mg/L)	Concentration (mg/L)	Concentration (mg/L)
Glycine	10.0	10.0	10.0	10.0
L-Arginine	200.0	200.0	200.0	200.0
L-Asparagine	50.0	50.0	50.0	-
L-Asparagine•H ₂ O	-	-	-	56.82
L-Aspartic acid	20.0	20.0	20.0	20.0
L-Cystine • 2HCl	65.0	65.0	65.0	65.2
L-Glutamic Acid	20.0	20.0	20.0	20.0
L-Histidine	15.0	15.0	15.0	15.0
L-Hydroxyproline	20.0	20.0	20.0	20.0
L-Isoleucine	50.0	50.0	50.0	50.0
L-Leucine	50.0	50.0	50.0	50.0
L-Lysine hydrochloride	40.0	40.0	40.0	40.0
L-Methionine	15.0	15.0	15.0	15.0
L-Phenylalanine	15.0	15.0	15.0	15.0
L-Proline	20.0	20.0	20.0	20.0
L-Serine	30.0	30.0	30.0	30.0
L-Threonine	20.0	20.0	20.0	20.0
L-Tryptophan	5.0	5.0	5.0	5.0
L-Tyrosine disodium salt dihydrate	29.0	29.0	29.0	28.83
L-Valine	20.0	20.0	20.0	20.0
Biotin	0.2	0.2	0.2	0.2
Choline chloride	3.0	3.0	3.0	3.0
D-Calcium pantothenate	0.25	0.25	0.25	0.25
Folic Acid	1.0	1.0	1.0	1.0
Niacinamide	1.0	1.0	1.0	1.0
Para-Aminobenzoic Acid	1.0	1.0	1.0	1.0
Pyridoxine hydrochloride	1.0	1.0	1.0	1.0
Riboflavin	0.2	0.2	0.2	0.2
Thiamine hydrochloride	1.0	1.0	1.0	1.0
Vitamin B ₁₂	0.005	0.005	0.005	0.005
i-Inositol	35.0	35.0	35.0	35.0
Calcium nitrate (Ca(NO ₃) ₂ • 4H ₂ O)	100.0	100.0	100.0	100.0
MgSO ₄ (anhyd.)	48.84	48.84	-	48.80
Magnesium Sulfate (MgSO ₄ • 7H ₂ O)	-	-	48.84	-
Potassium Chloride (KCl)	400.0	400.0	400.0	400.0
Sodium Chloride (NaCl)	6000.0	6000.0	6000.0	6000.0
Na ₂ HPO ₄ anhydrous	800.0	800.0	-	800.7
Sodium Phosphate dibasic	-	-	800.0	-
Sodium Bicarbonate (NaHCO ₃)	1500.0	2000.0	2000.0	2000.0
D-Glucose (Dextrose)	4500.0	2000.0	2000.0	2000.0
Glutathione (reduced)	1.0	1.0	1.0	1.0
L-Glutamine	300.0	300.0	300.0	300.0
Phenol Red	5.0	5.0	5.0	-
HEPES	2383.0	-	5958.0	-
Sodium Pyruvate	110.0	-	-	-
Penicillin G Sodium salt	-	-	67.0	-
Streptomycin Sulfate	-	-	100.0	-